

Ilastik classification, tracking, and manual editing workflow

Braden Stephens Keiser

2022/06/13



Pixel Classification

Tracking with Learning (Pixel Prediction Map)

Manual Tracking Editing with Custom Python Script

Install Ilastik from this file folder

Install TrackMate CSV Importer




Puncta Pixel Classification

Create New Project

- Pixel Classification
- Autocontext (2-stage)
- Pixel Classification + Object Classification
- Object Classification [Inputs: Raw Data, Pixel Prediction Map]
- Object Classification [Inputs: Raw Data, Segmentation]
- Manual Tracking Workflow [Inputs: Raw Data, Pixel Prediction Map]
- Tracking [Inputs: Raw Data, Binary Image]
- Tracking [Inputs: Raw Data, Pixel Prediction Map]
- Animal Tracking [Inputs: Raw Data, Binary Image]
- Animal Tracking [Inputs: Raw Data, Pixel Prediction Map]
- Tracking with Learning [Inputs: Raw Data, Binary Image]
- Tracking with Learning [Inputs: Raw Data, Pixel Prediction Map]
- Carving
- Boundary-based Segmentation with Multicut
- Cell Density Counting
- Data Conversion

Open Project ...







 Browse Files

Start a new project and specify directory

Import your image file

Change Z to T by 're-order hyperstack' in ImageJ

Then, your image should be 'tyx' (if 1 channel)




ERROR (in tracking) if you do not change Z→T

Choose **features**:


This relates to the size of your object


Because the **puncta are small**, we choose about just the first half of sigma values: **0.30 – 1.60**

If using **cells**, use **3.50-10.0**





Examine the selected features here. This is what the computer will be using for learning.









Change color of your objects by double-clicking the color
Paint the area that is your **specific object**
Paint the **background** **everywhere else**




Click 'live update' to see the first learning











Live Update:
You can keep this always on.



Switch between original
and annotated: i





Export:
 'Probabilities' = 'Pixel-Prediction Map'
 'Simple Segmentation' = 'Binary Image'

Hdf5 = BEST format for ilastik (and large datasets)

Tiff sequence = only good for viewing image in imageJ—exports the image ONE by ONE

Name	Status	Date modified	Type
base atg8 with t_Probabilities_00	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_01	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_02	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_03	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_04	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_05	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_06	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_07	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_08	✓	6/13/2022 10:02 AM	TIFF File
base atg8 with t_Probabilities_09	✓	6/13/2022 10:02 AM	TIFF File

HDF5 export format

HDF5 is a flexible cross platform binary data format. Using Python, you can access the data inside the **h5** files using the **h5py** library.

The **h5** file contains two items at the root level, **images** and **table**. Using **h5py**, you can access them like a python dictionary:

```
import h5py
data = h5py.File("path/to/data.h5", "r")
data.keys()
# <KeysViewHDF5 ["images", "table"]>
```

In the **images** group you will find a subgroup for each object (identifiable by **object_id**) which contains two datasets: **labeling** (cutout of the segmentation showing the object) and **raw** (cutout of the raw data showing the object).

```
data["images"].keys()
# <KeysViewHDF5 ["0", "1", "2", "3", ...]>
```

Access individual the **raw** and **labeling** images for each object using the key for the object:

```
data["images"]["1"].keys()
# <KeysViewHDF5 ["labeling", "raw"]>
data["images"]["1"]["labeling"]
# <HDF5 dataset "labeling": shape (28, 26, 13), type "|u1">
# stored as a numpy array, access the data using numpy indexing:
data["images"]["1"]["labeling"][:1
```

[ilastik - Object Classification](#)

<https://www.ilastik.org/documentation/objects/objects.html>

(ctrl + f : 'hdf5')

This is not necessary to use now.



Later, we can see if it's helpful.


Working with CSV is easiest and convenient for all lab members.

Tracking with Learning

Import your previous image
Import your exported
prediction/probability image

- Manual Tracking Workflow [Inputs: Raw Data, Pixel Prediction Map]
- Tracking [Inputs: Raw Data, Binary Image]
- Tracking [Inputs: Raw Data, Pixel Prediction Map]
- Animal Tracking [Inputs: Raw Data, Binary Image]
- Animal Tracking [Inputs: Raw Data, Pixel Prediction Map]
- Tracking with Learning [Inputs: Raw Data, Binary Image]
- Tracking with Learning [Inputs: Raw Data, Pixel Prediction Map]
- Carving
- Boundary-based Segmentation with Multicut
- Cell Density Counting
- Data Conversion







You will see an “input” selection. These are your background and puncta. We will only focus on the puncta color.

I started a new test, so my atg8 is now yellow instead of red.



Select a **low minimum value**. These **puncta** are only a few pixels. So, I set the min to just **3 pixels**. Then, click “apply”, and we are finished with this part.


SKIP “Division Detection”
GO TO “Object Count Classification”





Select a few objects, and then do the “live update”.
Now, we can move to “training”.







Label your tracks frame by frame
It is good to make about **10-20 tracks** to increase the accuracy of predictions

In 'Tracking', Ilastik will not use these like a binary image. The program will make new tracks based on these training data.

TIP:

press 'S' to start a new track
Press 'E' to turn off labeled tracks
Press 'R' to turn off prediction
'I' has no function here



6. Tracking - Structured Learning

Detection	0.60000
Division	0.60000
Transition	0.01000
Appearance	0.30000
Disappearance	0.20000

Calculate Tracking Weights

Track!

Solver: Flow-based

Tracking Parameters:

Transition Neighbourhood	1
Max. Distance	200
Max. Object per Merger	1
Timeout (sec.)	

7. Tracking Result Export

1st: Turn **OFF**
'Divisible
Objects' or
the tracking
won't work

Tracking Parameters:

Transition Neighbourhood	1
Max. Distance	200
Max. Object per Merger	1
Timeout (sec.)	
Border Width (abs.)	10
<input checked="" type="checkbox"/> Divisible Objects	

UNCHECK!

Tracking Parameters:

Transition Neighbourhood	1
Max. Distance	200
Max. Object per Merger	1
Timeout (sec.)	
Border Width (abs.)	10
<input type="checkbox"/> Divisible Objects	

Detection	0.60000
Division	0.30000
Transition	0.80000
Appearance	0.60000
Disappearance	0.50000

Calculate Tracking Weights

Track!

Solver: ILP

Tracking Parameters:

Transition Neighbourhood	1
Max. Distance	25
Max. Object per Merger	1
Timeout (sec.)	

These 4 parameters have the most effect.

Appearance: if this is large, it will be harder to make new tracks

Disappearance: if this is large, it will be harder to end a track (it will try to 'survive')

Solver: The type of algorithm for calculating the tracks. This worked better for me than 'Flow' with puncta.

Max. Distance: This is the max distance that the puncta can travel between frames before it will be seen as a new track.

Example:

The cell size is about 30-40 pixels

The max. distance is **at most 40**

This is usually smaller when we look at the activity in each cell

This effect seems less than **above**
But: the **larger** of these **three values** will have the **strongest** effect

Detection: > Division/Transition will favor detection

Transition: > detection/division will favor transition (change current track to new track)


Division: we don't use it because no cells are dividing



e.g.: **transition > detection > division**

*****Transition** helps separate nearby tracks—adjusting this can be helpful



BEFORE EXPORT:


Open original image (**with T**) and export with BigDataViewer
as XML/HDF5—choose the location. **File = export.xml**






Export Tracking Results:
Source: PLUGIN
Plugin: Fiji-MaMuT
BDV file: export.xml of the original image





base atg8 with t_Fiji-MaMuT_mamut



TrackMate on export2

Display options

Edit settings

☒ Display spots ☒ as ROIs

Spot display radius ratio:

Display spot names: ☐

Color spots by:

auto min max

☒ Display tracks

Fade tracks in time: ☒

Fade range: time-points

Color tracks by:

auto min max

☐ Limit drawing Z depth pixel

TrackScheme **Tracks** Spots

Save Next

Track tables

Export to CSV coloring

	Label	Spot ID	Track ID	Quality (quality)	X (pixel)	Y (pixel)
Spots						
Edges						
Tracks	track-16	91	16	3	108.2	190.6
	track-16	93	16	3	86	128.2
	track-16	14	16	3	173.4	205.6
	track-16	62	16	3	119.148	193.296
	track-16	78	16	3	118.568	199.946
	track-17	97	17	3	218.444	149.333
	track-17	37	17	3	270.615	187.615
	track-17	102	17	3	217.467	149.333
	track-17	71	17	3	240	248.167
	track-17	109	17	3	273.125	186.875
	track-17	15	17	3	242.385	231.846
	track-17	79	17	3	240.333	248.667
	track-17	111	17	3	118	82
	track-17	48	17	3	270.077	188.615
	track-17	113	17	3	124.667	81.75

```

dir_in = 'C:/Users/Braden/OneDrive/Documents/Codings/machine Learning/Lihui atg2 atg8/testing week 3 0530/coding'
path_atg8 = os.path.join(dir_in, 'export atg8 0606.csv')
atg8 = pd.read_csv(path_atg8)

dis = 11

df = atg8[['LABEL', 'FRAME', 'POSITION_X', 'POSITION_Y']]
df = df.drop([0,1,2], axis = 0)

df['POSITION_Y'] = df['POSITION_Y'].astype(float).astype(int)
df['POSITION_X'] = df['POSITION_X'].astype(float).astype(int)
df['FRAME'] = pd.to_numeric(df['FRAME'], downcast = 'unsigned')

df['puncta']=None

def find_jumpers(df):

```

lihui atg2 atg8 > testing week 3 0530 > coding

Name	Status	Date modified
part 1 atg8 presentation		6/6/2022 4:41 PM
part 1 original and jumpers functions.py	✓	6/6/2022 4:41 PM
export atg8 0606	✓	6/6/2022 4:38 PM
part 1 original and jumpers functions ilastik ...	✓	6/6/2022 11:27 AM
part 2 colocalization function.py	✓	6/6/2022 11:06 AM

Search coding

Copy address
Copy address as text
Edit address
Delete history

Max. distance track can travel in the next frame

```

for x in df:
    find_jumpers(df, 23)
    if 'maybe' in df['puncta']:
        find_jumpers(df, 23)
#df.to_csv(os.path.join(dir_in, 'part 1 atg8 presentation 0606.csv'), index = False)

```

```

return invalid
for x in df:
    find_jumpers(df, 23)
    if 'maybe' in df['puncta']:
        find_jumpers(df, 23)
df.to_csv(os.path.join(dir_in, 'part 1 atg8 presentation 0606.csv'), index = False)

```

Name	Status
part 1 atg8 presentation 0613	✓
part 1 original and jumpers functions addin...	✓
part 1 original and jumpers functions.py	✓
part 1 atg8 presentation 0606	✓
export atg8 0606	✓

base atg8 with.tif

1/20, 275x275 pixels; 8-bit, 1.4 MB


Color Inspector 3D
DeconvolutionLab2
Examples
Feature Extraction
HDF5
Image5D
Integral Image Filters
Janelia H265 Reader
LOCI
LSM Toolbox
Landmarks
MaMuT
Multiview Reconstruction
Optic Flow
PSF Generator
Process
Registration
Ridge Detection
SPIM Registration
Segmentation
Skeleton
Stitching
Time Lapse
Tracking

MTrack2
Manual Tracking
ToAST
TrackMate
Load a TrackMate file
Manual tracking with TrackMate
TrackMate CSV importer
FilamentDetector


23
LABEL
FRAME
POSITION_X
POSITION_Y
puncta
TRACK_ID
POSITION_T

In [30]: run
0530/coding/0
23


atg8 presentation 0606.csv', index = False)



Make sure this is an integer with no decimal value!




Original, no function




Probably max.
distance
problem


With function
(original tracking parameters)



Original, no function



With function
(with modification to tracking
parameters)



With function
(with modification to tracking
parameters)

Index	LABEL	FRAME	POSITION X	POSITION Y	puncta	TRACK ID	POSITION T
3	track-32	18	145	33	original	32	18
4	track-32	19	134	31	original	32	19
5	track-33	18	71	114	original	33	18
6	track-33	19	73	119	original	33	19
7	track-2	0	31	1	original	2	0
8	track-43	11	122	93	from track-43	43	11
9	track-39	6	137	52	from track-39	39	6
10	track-36	12	139	91	from track-36	36	12
11	track-39	7	135	53	from track-39	39	7
12	track-48	3	130	38	from track-48	48	3
13	track-36	13	139	92	from track-36	36	13

```

data['POSITION_T'] = data['FRAME']
print(tracked)

find_original(df)

def find_jumpers(data, tracked):
    jumpers = data[data['puncta'] != 'original']
    label = jumpers.groupby(by = ['LABEL'])
    size = tracked + 10
    for tracks, series in label:
        track = series.sort_values(by = 'FRAME', ascending = True)
        x, y, t = 'POSITION_X', 'POSITION_Y', 'FRAME'
        indice = track.loc[track['FRAME'] == track[t].min()]
        origy = int(track[y][indice.index])
        origx = int(track[x][indice.index])
        x_range = range(origx - dis, origx + dis, 1)
        y_range = range(origy - dis, origy + dis, 1)
        names = tracks
        for index, row in track.iterrows():
            if row[x] in x_range and row[y] in y_range:
                data.at[index, 'puncta'] = f'from {names}';
                data.at[index, 'LABEL'] = f'track-{size}'
            elif row[x] not in x_range and row[y] not in y_range:
                data.at[index, 'puncta'] = 'maybe'
            elif row[x] not in x_range or row[y] not in y_range:
                data.at[index, 'puncta'] = 'maybe'
            else:
                return 'invalid'

```

The dataframe created with the functions fulfills the TrackMate requirements.

'Puncta' is helpful to the user. Original = the tracks Ilastik originally found.

'From--##' tells us if ilastik was wrong. If ## < 'size', then this is the first track that the function gives to new puncta after separating the tracking error from the original. If ## > 'size', then this was a second puncta tracking error, and the function created a *second new track*.

Other CSV filetypes:

```

FOR OTHER CSV FILES:
x = 'POSITION_X'
y = 'POSITION_Y'
t = 'FRAME'
FRAME = 'FRAME'
TRACK_ID = 'LABEL'
name = 'LABEL'

and then you can perform the following:
df = pd.DataFrame(columns = [x, y, t, name, TRACK_ID, FRAME, 'puncta', 'LABEL'])
df[x] = atg2['CENTER_OF_OBJECT_0'] etc...

FOR POSITION_X, POSITION_Y, POSITION_T, FRAME we need to CONVERT from STRING to INT:
workflow:
df[x] = atg2['CENTER_OF_OBJECT_0']
df[x] = df[x].astype(float).astype(int)
etc...

you then transforms the column data from the original CSV into a format which is more readable to trackmate
i.e. track_in , position_t, position_x, etc...
trackmate needs:
track_id
position_x
position_y
position_t (ideally) or z
frame

adding 'puncta' and 'label' is useful for checking the progress while in the python IDE

```

You can use any tracking file; however, you just need to specify the five basic parameter values that TrackMate needs:

POSITION_X
POSITION_Y
POSITION_T
FRAME
TRACK_ID

Ilastik CSV export (instead of MaMut)

Example:

x = 'POSITION_X'
y = 'POSITION_Y'
Etc...

```
df = pd.DataFrame(columns = [x, y, t, name, TRACK_ID, FRAME, 'puncta', 'LABEL'])
```

```
df[x] = input_data['Object_Center_0']
df[y] = input_data['Object_Center_1']
```

	N	O	
Object_Center_0	Object_Center_1	Ob	
31.608696	1.434783		
133.300003	47.933334		
144.399994	79.040001		
124	88.714287		
94.647057	132.529419		